



## Origami instructions for the BioBox™



Hold the base / lid with both hands and make sure it is smooth side down. Use your thumbs and palms to fold the two longest sides up at 90° from the back, while using your fingers to support the card.



Still holding the base / lid with the long sides folded up, bring the near side of the base / lid close to the body for support. Use the middle fingers to lift the short sides while using both index fingers to push in at the corners.



Bring the base / lid close to the body for support. Once the short sides have folded up bring the index fingers out of the corner folds to help support the sides in the up position.



Make sure the corner folds are folded flat against the sides. You can now use both thumbs to fold the rest of that side down into the box making sure the inner corners are pushed into lock position.



Turn the base / lid 180° and repeat all the steps above.